

# Data acquisition developments and collaborations at Daresbury Laboratory

# Collaboration

- Collaboration not necessarily easy
- Communication is key
- Best to practise with friends first
- Will describe couple of developments based on collaborations
  - PXGEN++
  - DNA

# PX Objectives

- PXGEN out of date and hard to update.
- Robust replacement offering greater functionality.
- Uniform interface across stations, detectors, cameras etc.
- Simple GUIs that make set up and data collection easy.
- Automation wherever possible.
- Link data acquisition and processing.
- Flexibility to cope with future changes.
- Remote monitoring and control.

# Daresbury Objectives

- Cross facility group collaboration.
- Possible deployment to DIAMOND.
- Plug and play construction.
- Scripting language to facilitate easy configuration of experiments.
- Distributed solution.

# Strategy


- Inter project collaboration.
  - PX, NCD, XRS, XRD, VUV-IR.
- Incremental and iterative development.
- Largely object oriented development. Java, C++ and C.
- Classes to represent real beam line hardware e.g. detectors, cameras, tables, mirrors etc.
- Use interfaces to give plug and play.
- Use CORBA to enable distribution.
- Provide common ‘look and feel’ where appropriate.
- Configuration of objects via XML files.

# Advantages to in-house collaboration

- Shared development effort allows more to be achieved without repetition.
- Other programmers understand the software.
- Wider use leads to greater robustness.
- Get benefits from aspects that may not at first be seen as directly PX aims e.g. scripting, XML configuration.
- Users see common interfaces across beam lines and stations.
- Aids the development of multiple technique stations.


# PXGEN++

- Same GUI for data collection independent of detector in use.
- Data collection ready, tested with Q4 and MAR image plates (MAR CCD next).
- Using OEMove for beam line set-up.
- PXGEN++ and OEMove undergoing final tests.
- Once fully tested and released on one PX station should be able to deploy on others.
- Need to integrate in XAFS control, automation, sample changing etc.
- Standard XAFS GUI and control under development as part of 6.2 project (NCD, XRD and XAFS).





# Detectors and Cameras



Abstract Detector  
provides standard  
interface e.g.


- initialise()
- clear()
- read().....

Concrete detectors  
implement these.

Similarly for Camera.

GUI/Experiment  
only know that they  
have a Detector and  
a Camera.

# OEMove



- OEMove provides a standard tool to allow the control and set up of Optical Elements.


- OEedit allows the configuration of one or more visual representations of OEs.

- OEMove reads in representations at run time.

- One tool meets needs of PX and other stations.

# Optical Elements


- OEs obey a standard interface.
- Each OE has a number of specific DOFs (Degree Of Freedom).
- Specific DOFs are implemented generally e.g. FrontHoriz and BackHoriz by SingleAxisLinearDOF.
- DOFs obey a standard interface.



# Scripting - JCLAM

- JCLAM is a scripting language developed on top of JPython. A set of Python scripts basically.
- Can still use standard Python language.
- JCLAM run separately and accepts commands via a message handler.
- Separates control from the GUI.
- Standard scanning GUI will be set up.
- Allows easy/rapid alteration of control.
- Users can be given control of hardware and allowed to alter scans where required e.g. in different projects.
- Examples
  - pos LeftSlit 2.0 RightSlit 2.0
  - scan Slit 1 1 5 qcr 1

# Absorption Edge Scans



- Standard absorption edge scan interface.
- Supports multiple modes, Angle, KSpace and Fast.
- GUI completely detached from control.
- Scans done via JCLAM.

`scansie.SCANSIE(['mono',18235.3,-1.8,18054.2,'tm',1000.0])`



# Automation

- Development of Optical Elements and the use of JCLAM should facilitate beam line automation.
  - Feedback required on beam
  - Objectives include
    - Alignment of camera to beam, standard peak searching and scanning.
    - Slit optimisation.
    - Focussing of mirror/monochromators.
    - Sample changing and alignment.
- Feedback from data processing (DNA).
- Automatic logging of data collection.

# DNA - DNA's Not Autostruct

- A collaboration presently involving Cambridge LMB, ESRF and Daresbury.
- Still at a relatively early stage and under active development (unfortunately part time).
- Eventual aim is to completely automate the collection and processing of Protein Crystallography data.
- Initially the project will develop what we are calling an 'expert system' sitting between data collection software at a given x-ray synchrotron and Mosflm.
- Utilising the Mosflm server under development at Cambridge for new Mosflm GUI.
- The project however also aims to allow the use of other data processing software in the long term.
- <http://www.dna.ac.uk/>

# DNA architecture





# DNA progress at ESRF

- ProDC modified to add a button, characterise crystal.
- Expert system developed to communicate with Mosflm server.
- ProDC communicates with Expert system.
- Communication to/from expert system and Mosflm server work.
- Crystal mounted on ID14.
- Characterisation collected two images on behalf of expert system.
- Expert system instigated auto-indexing of images by Mosflm via Mosflm server.
- Results relayed back to ProDC and displayed in message window.

ProDC Beamline Control ProDC V7.6b

File View Commands Database Help

Instrument(s) Control

Fast X-ray shutter is **CLOSED**

Motor Position Control:	Current	Destination	Relative
Xtal-Det dist	1.35	0.0	0.0
Phi	+90 -90 6.0	0.0	0.0
Mono Bragg	5.0	0.0	0.0
Wavelength	0.0	0.0	0.0
Energy	0.0	0.0	0.0

Darren's Detector

File format:

Save Corrected Image

No Image Display

☐ Take dark current at start of collection

Binned Fast

Status Messages

Data Collection Input

Collection mode **Single Run**

Exposure Mode **Oscillation**

Frame number to start	2	Directory	/data/id14eh2/inhouse/dna/
Oscillation Start (deg)	90.000000	Prefix	ntest3
Oscillation Range (deg)	1.000000	Run number	1
Overlap between frames (deg)	0	Comments	1.47
Exposure time per pass (sec)	10.000000	Template	ntest3_1_###.img
No. of passes per frame	1	Disk space (kb)	1000
No. of frames to collect	1	Total exposure time (sec)	10.0
		Total execution time	13.7 seconds

Start data collection

Stop data collection

Characterise Crystal

Enable Button Box

Abort Now

Detector Distance (m): 1.35  
 Wavelength: 0.933  
 Data Directory: /data/id14eh2/inhouse/dna/  
 Filename suffix: img

09:04:25 AM: Frame 2 : 90.0 to 91.0 degrees

Collection finished\_09:04:25 AM

Result of auto index:

Space group: F23

Cell: a= 180.76, b= 180.76, c= 180.76, alpha= 90.00, beta= 90.00, gamma= 90.00

Updated x\_beam= 93.61, y\_beam= 91.98

Reflection used=1331, rejected= 53

Spot deviation (RMS) 0.100

# DNA progress at DL

- Stand alone GUI Built
- Can use GUI from PXGEN++
  - Will contain results of strategy calculation and allow modification before data collection
- Expert system used as developed at ESRF
- Java class written to interface to expert system
- Communication to/from same expert system works
- XML parsing via castor, <http://castor.exolab.org>  

```
Object o = unmarshaller.unmarshal(new StringReader(command));  
CollectStart cs = (CollectStart) o;
```

**PXGEN++ (V1.1)**

**Quit**

**Setup**

Beam Alignment    Beamline    Goniometer    **Detector Distance**

Detector Specific    Fluorescence Scan



**Collect**

Edit Exposures    Perform Exposures    Single Exposures

**Abort Exposures**    Expert System

**Status**

**Shutter**

**Phi**

**91.0**


**Sample to Detector**

**250.0**

**Operation**

ready

**Progress to End**



Mar 13, 2002 14:04

**Expert System**

**Data directory**    /usr/users/shk/mosflm-server/server/data

**Prefix filename**    test\_

**Detector Distance**    250    mm

**Wavelength**    0.979    Angstrom

**Number of Experiments**    2    ▼

Start Phi	Phi Increment	Number of Images	Time	Template Number
10.0	1.0	2	5.0	0
30.0	2.0	5	2.0	10

Image ready, here you go clever clogs...

Result of auto index:  
 Space group: P4  
 Cell: a= 79.19, b= 79.19, c= 37.92, alpha= 90.00, beta= 90.00, gamma= 90.00  
 Updated x\_beam= 94.01, y\_beam= 91.88  
 Reflection used=1100, rejected= 88  
 Spot deviation (RMS) 0.224

**Characterise Crystal**    **Collect**    **Dismiss**

# Acknowledgements

## Daresbury

Colin Nave, Liz Duke, Geoff Mant, Paul Stephenson, Mike Miller, Karen Ackroyd, Dave Love

## ESRF

Olof Svensson, Darren Spruce

## Cambridge

Andrew Leslie, Harry Powell, Graeme Winter